

VISOCOLOR® ECO

Testanleitungen

Test instructions

Instructions des tests

Instrucciones de los tests



zur visuellen Bestimmung
for visual determination
pour la détermination visuelle
para la determinación visual


visicolor® ECO Gesamthärte

total Hardness / Dureté totale / Durezza total


Messbereich / Range / Domaine de mesure / Rango

1 °d \triangleq 1 °d


REF 931 029




1
5 mL Probe
Sample
Echantillon
Muestra




2
2 °d GH-1



3
Umschwenken
Shake gently
Secouer légèrement
Mezclar volteándolo



4
°d GH-2
bis Farbumschlag
till coloration
jusque inversion de couleur
hasta cambio de color



5
1 °d \triangleq 1 °d
rot → grün
red → green
rouge → vert
rojo → verde



1:30

Meerwasser
Sea water
Eau de mer
Agua de mar